

Handheld 2

Truly portable tube reader

The Ziath Handheld 2 reader is a truly portable reader which can be held and operated in one hand; leaving the other free for tube handling.

Features

- Compatible with all tubes
- Portable, independent and compact
- No external computer or tablet needed

Easy to Set-Up and Easy to Use

The Handheld 2 tube reader requires no external computer and is easily operated with one hand; leaving the other hand free to manipulate tubes. The software is pre-loaded and boots up immediately with a battery life of up to a week with casual use. This device allows you to read tubes standing in the basement at the freezer or liquid nitrogen tank. Even out in the field acquiring samples.



Handheld 2

Truly portable tube reader

Operation Modes

The reader has four operation modes:

- **Barcode reading** - Read a tube barcode to view it on the screen
- **Picklist** - Enter a picklist using the associated Client and read the tubes with the handheld device to ensure you have the correct tubes while at the freezer
- **Data View** - Read a barcode to see the details of the sample in the tube
- **Data Entry** - Read a barcode and enter the information on the same to store

Client App

To assist with the transfer of data on and off the handheld so that it can be used in environments with no computer network available, the handheld comes with an associated client application to help define the data to store and load/unload it onto the handheld. This allows you to go into the deepest darkest corners of your sample storage and still be able to access your data!



Ordering information

Code	Description
601090	Handheld 2 Single Tube 2D Barcode Reader



Ziath Ltd Email: info@ziath.com www.ziath.com
An Azena Life Sciences company www.azenta.com

Ziath technology is protected by patents and patent applications in the UK, USA and worldwide. Ziath has a continuous development programme aimed at further improving and developing its products. All specifications are, therefore, subject to change. E&OE. All trade marks acknowledged. © 2023 Ziath Limited. 44012-SO 0923



Biozym Scientific GmbH
Tel.: 05152 / 9020, Fax: 05152 / 2070
Mail: support@biozym.com

Biozym Biotech Trading GmbH
Tel.: 01 / 334 0156 0, Fax: 01 / 334 0156 88
Mail: support@biozym.com



Biozym Scientific



biozym_scientific



Biozym Scientific



Newsletter